

Display Report

Analysis Info

Analysis Name D:\Data\test\0401\moshizisuan_4_01_3137.d
Method 3137.m
Sample Name moshizisuan
Comment

Operator
Instrument brucker
amaZon SL

Acquisition Parameter

Ion Source Type	ESI	Ion Polarity	Positive	Alternating Ion Polarity	on
Mass Range Mode	UltraScan	Scan Begin	70 m/z	Scan End	1000 m/z
Capillary Exit	140.0 V	n/a	n/a	Trap Drive	50.6
Accumulation Time	1892 μ s	Averages	5 Spectra	Auto MS/MS	off

